
Cefar URO

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TABLES DES MATIÈRES

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
User Manual

CONTENTS

1. Introduction	50
2. How to use the medical equipment (Intended Use)	51
2.1 Fields of Application	51
2.2 Therapy objectives	51
2.3 Contraindications	52
3. Warnings	53
4. Overview	59
4.1 Device components and accessories	59
4.2 Explanation of symbols	59
5. Description of the Cefar URO	61
5.1 Description of the device	61
5.2 Explanation of symbols	62
5.3 Battery Charge	63
5.4 Accessories	64
6. Step by step use - Patient mode	65
6.1 Attach the electrodes	65
6.2 Switch on the stimulator	65
6.3 Select a preset Program	66
6.4 Adjust timer settings if needed	66
6.5 Start the stimulation	67
6.6 Key Lock	67
6.7 Program pause	67
6.8 Open lead	68
6.9 Resume treatment after open lead	68
6.10 Program end	69
7. Physician Mode	70
7.1 To access Physician Mode	70
7.2 Physician Mode home screen	70
7.3 Lock Program	71
7.4 Custom Programs	71
7.4.1 Select custom Programs	71
7.4.2 Custom Program first screen	71
7.4.3 Summary screen	72
7.4.4 Parameters range	72

7.5 Compliance.....	73
7.6 CC-CV.....	73
7.7 Pain score.....	73
8. Programs	74
8.1 Program parameters.....	74
8.2 Program description and recommendations.....	75
9. Troubleshootings.....	76
9.1 Electrode or cable failure.....	76
9.2 Battery level.....	76
9.3 Others.....	77
10. Care, Maintenance, Transport, Environmental Statement	79
10.1 Maintenance.....	80
10.2 Transport.....	80
10.3 Environmental Statement, Expected Life	81
11. Technical Datas, Standards, Guarantee, Patients	82
11.1 Technical Data.....	82
11.1.1 General information	82
11.1.2 Stimulation Parameters.....	82
11.1.3 Information on electromagnetic compatibility (EMC)	83
11.1.4 Environmental conditions.....	83
11.2 Standards.....	84
11.3 Guarantee	84
12. Frequently Asked Questions	85
13. Contacts.....	86
14. Electrode placement guide.....	87

1. INTRODUCTION

As part of Chattanooga®, CEFAR stimulators are widely used and endorsed by medical professionals around the world.

TENS is effective and economical. Through clinical research, areas of application for Transcutaneous Electrical Neurostimulation (TENS) and Percutaneous Neurostimulation (PENS) are rapidly expanding. Chattanooga®/Cefar is working actively to further develop the method towards a natural treatment alternative for the patient.

Cefar URO is a dual channel nerve stimulator (TENS) designed for use by both professionals and home users. It inhibits detrusor contraction by modifying the integration of afferent and efferent information. The stimulator features 4 preset programs and 3 custom programs. Program and amplitude can be individually set for each channel, which means that two different programs can be combined simultaneously.

Chattanooga® is the world's largest manufacturer of rehabilitation and pain relief equipment. Since 1947,, Chattanooga® products have contributed to treating patients in hospitals, clinics and home settings worldwide.

Chattanooga® is a brand of ENOVIS™. ENOVIS™ provides solutions for musculoskeletal health, vascular health and pain management. Our products help prevent injuries or rehabilitate after surgery, injury or from degenerative disease and relieve pain enabling patients to regain or maintain their natural motion. Visit enovis-medtech.eu

2. HOW TO USE THE MEDICAL EQUIPMENT (INTENDED USE)

Note

- This manual is considered as an accessory of the therapy unit and therefore it should accompany it at all times.
- The specific instructions provided here are conditions for the intended use and correct operation of the equipment as well as the safety of the patient and the operator using it.

2.1 Fields of Application

The Cefar URO is a stimulator used by patients at home to ensure electric stimulation treatments to regulate bladder instability.

The Cefar URO is a stimulator used by medical professionals to apply neuro muscular stimulation (EMS/ NMES) treatments

TENS = Transcutaneous Electrical Nerve Stimulation. TENS is an effective method offering treatment of bladder instability. TENS stimulates the Posterior Tibial Nerve to generate neuromodulation at the sacral level, thus treating bladder instability.

2.2 Therapy objectives

As a TENS device, indications are for treatment of bladder instability (overactive bladder).

With the Pain Relief programs, impulses are sent through the nerves to generate a neuromodulation.

The Cefar URO is a electrotherapy unit that can be used to treat bladder instability. All the programs can be used as often as required.

Guidelines for using the preset treatment programs with indication examples, suggested electrode placements and recommended treatment time are included in this manual.

2.3 Contraindications

Danger!

Do NOT use the Cefar URO on patients with:

- Implanted electronic devices. Do not use the device if you have a cardiac stimulator, implanted defibrillator or other implanted electronic/electrical device. This could cause an electrical shock, burns, electrical interference or death.
- Epilepsy
- Pregnancy (do not use on abdominal region)
- Serious arterial circulation problems in lower limbs
- Do not use chest stimulation on patients with cardiac arrhythmia or angina pectoris
- Heart disease: If you have suspected or diagnosed cardiopathy you should follow the precautions for use recommended by your doctor because the device may cause lethal rhythm disturbance to the heart in susceptible individuals.

Note

Fixation metalwork and/or prosthesis equipment

The presence of fixation metalwork and/or prosthesis equipment (metallic equipment in contact with the bone: pins, screws, plates, prostheses, etc) is not a contraindication for TENS stimulation. The electrical current of the Cefar URO is specially designed to have no harmful effect on osteosynthesis equipment.

Electro-stimulation should not be applied if any part of the implanted element going through the patient's skin because this may cause shocks or burns.

3. WARNINGS

Definitions

It is mandatory to read the safety statements before using the physiotherapy unit. The safety statements are classified as follows:

**Danger!**

This term indicates an imminent hazard. If not avoided, this hazard could result in death or serious injury.

**Warning!**

This term indicates a hazard. If not avoided, this hazard can result in death or serious injury.

**Caution!**

This term indicates a potential hazard. If not avoided, this hazard can result in minor personal injury and/or product/property damage.

Safety Information

**Danger!**

Explosion hazard - Cefar URO is not designed for use in areas where an explosion hazard may occur. An explosion hazard may result from the use of flammable anaesthetics, oxygen-rich environments, skin cleansing agents and disinfectants.

**Warning!**

Patient hazard

- Only authorised individuals are allowed to repair the Cefar URO. Individuals are authorised after receiving training in the operation of the unit and reading this operating into the manual.
- Before using the therapy unit, the operator must ascertain that it is in correct working order and operating condition. The cables and connectors, in particular, must be checked for signs of damage. The device cannot be used if any part is damaged as there is a risk of electrical shock.

- Stop therapy immediately if you have doubts about the device settings and/or the therapy protocol.
- Patients must be fully conscious while being instructed in the use of the therapy unit and during therapy.
- The choice of the therapy parameters to program and of the therapy protocols to use is restricted to the responsible physician or therapist. It is the physician's or therapist's decision whether or not to use the unit on a specific patient.
- The patients must be familiar with the functions of the Cefar URO device allowing them to stop therapy, if needed. Patients unable to operate the emergency stop function e.g. paralytic patients, must never be left unattended during therapy.
- Any accessories used with Cefar URO must first be approved by the manufacturer. Using of not specified accessories and detachable parts could be unsafe for the user.
- The utmost caution is advised under the following conditions. Depending on the judgement of the responsible physician, the unit may only be applied under supervision and with the parameters defined by the responsible physician. Otherwise the exercise may be too strenuous for the patients with:
 1. Hypertension (> stage 2), ischemic heart disease and cerebrovascular diseases
 2. Cardiovascular diseases
 3. Pregnancy

Never apply the electrodes:

- Near the head or on the eyes
- On the front and side of the neck (especially carotid sinus) because this could cause severe muscle spasms resulting in closure of the airway, difficulty in breathing, or adverse effects on heart rhythm or blood pressure.
- Counter-laterally, i.e. do not use two poles connected to the same channel on opposite sides of the body.
- On or near skin lesions of any kind (swelling, burns, irritation, eczema, cancerous lesion, etc)
- Crossing over the heart because the introduction of electrical current into the chest might cause rhythm disturbance to the patient's heart, which could be lethal.

Warning!

- Keep the device out of reach of children
- Never leave the device unattended when it is switched on!
- After use, store the device in a safe place to avoid other people not informed to use the device.
- This device is not a toy but a medical device. Misunderstanding its use can cause damage!

Warning!

Shock hazard - Strictly observe the following warnings. Failure to do so could endanger the lives of the patient, the user and other persons involved.

- **Before use** allow the Cefar URO to reach room temperature. If the unit has been transported at temperatures below 0 °C (32°F), leave it to reach at room temperature for about 2 hours, until any condensation has disappeared.
- **Electrosurgical equipment or defibrillators.** Disconnect the electrodes from the device before using electrosurgical equipment, or a defibrillator, to avoid cutaneous burns from the electrodes and destroying the device.
- **Electronic surveillance equipment.** Do not apply stimulation near electronic surveillance equipment (e.g. cardiac monitors, ECG alarms), as there is a risk they may not work properly whilst the electrical stimulation device is being used.
- **Electromagnetic radiation.** Do not use the stimulator in areas in which unprotected devices are used to emit electromagnetic radiation. Portable communications equipment can interfere with the device.
- **Cancer.** Do not apply stimulation if you have progressive cancer or near any cancerous tumour. The increased metabolism, caused by certain modes of stimulation, is likely to encourage cancer cells to spread
- **Battery leakage.** If there is leak from a component, take steps to ensure the liquid does not come into contact with skin or eyes. Should this occur, wash the affected area with water and consult a doctor.
- **Strangulation.** Do not wind cables around the neck. Tangled cables can cause strangulation.
- **Post-surgery.** Proceed with caution after recent surgery when stimulation may disrupt the patient's healing process.
- **Accessibility of the power adaptor.** The plug socket must be close to the power adaptor and be easily accessible.
- The Cefar URO must only be operated in dry rooms.
- Do not use the Cefar URO in water or in a humid atmosphere (sauna, bath, shower etc.), etc.).
- When connecting the unit to other equipment or when creating a medical system, check that the sum of leakage currents will not cause any hazard. Please contact Enovis if you have questions regarding this matter.
- No modification of this equipment is allowed.
- Do not open the product and its accessories as there is risk of electrocution.
- Before cleaning and service interventions, **turn off the device.**
- **Liquids and foreign material (such as dust, metal etc) must not be allowed** to enter the device. If such material has entered into the units, it must be immediately checked by a service technician, before it can be reused.
- **Electricity supply.** Never connect the stimulation cables to an external power supply, as there is a risk of electrocution.
- Do not apply stimulation near the area of an implant, such as cochlear implants, pacemakers, skeletal anchorage or electric implants. This could cause an electrical shock, burns, electrical interference or death.

- Never use the Cefar URO or the AC adaptor if it is damaged or open. There is a risk of electric shock.
- Disconnect the AC adaptor immediately if there is abnormal heating or smell, or if smoke comes from the AC adaptor or the device.

Warning!

Equipment malfunction - this warnings can prevent equipment malfunctions that result in patient hazards

- Keep the unit away from such equipment and verify its performance before use.
- Do not use the Cefar URO within one metre of short wave or microwave devices as this could alter the currents generated by the stimulator. If you are in any doubt as to the use of the stimulator in close proximity to another medical device, seek advice from the manufacturer of the latter or from your doctor.
- Exercise caution when using electrotherapy while the patient is connected to monitoring equipment with electrodes attached to the body. Stimulation could disrupt the signals sent to the monitoring equipment.
- Refer repair and maintenance to authorized persons. Persons are authorized after training by a specialist trained and commissioned by the manufacturer.
- Inspect the Cefar URO and it's accessories for damage and loose connections at least once a year. Damaged and worn parts must be immediately replaced with original spare parts by authorized staff.
- Do not use the device if you are connected to a high-frequency surgical instrument as this could cause skin irritation or burns under the electrodes.

Caution!

Patient hazard - these cautions need to be observed to avoid the risk of electrical shock or other negative effects to the patient.

- Do not apply stimulation close to metal. Remove jewellery, piercings, belt buckles or any other metallic product or device in the area of stimulation
- Be careful if the patient has sensitivity problems or is not able to communicate that he or she feels discomfort, however light.

- Never begin an initial stimulation session on a person who is standing. The first five minutes of stimulation must always be performed on a person who is sitting or lying down. In rare instances, people of a nervous disposition may experience a vasovagal reaction. This is of psychological origin and is connected with a fear of the muscle stimulation as well as surprise at seeing one of their muscles contract without having intentionally contracted it themselves. A vasovagal reaction causes heart to slow down and blood pressure to drop, which produces a feeling weakness and a tendency towards fainting. If this does occur, all that is required is to stop the stimulation and for the person to lie down with the legs raised until the feeling of weakness disappears (5 to 10 minutes)
- Do not disconnect any channel during the stimulation.
- Do not use the stimulator while driving or operating machinery or during any activity in which electrical stimulation or involuntary muscle contraction can put the patient at risk of injury.
- Do not apply stimulation during sleep.
- Do not use the stimulator at altitudes of over 3,000 meters.
- When using constant current mode, always turn off the stimulator before moving or removing any electrodes during a session, to avoid electrical shock to the patient.
- Do not try to place electrodes on a body part not directly visible without assistance.
- Attach the electrodes in such a way that their entire surface is in contact with the skin.
- For obvious reasons of hygiene, each patient must have their own set of electrodes. Do not use the same electrodes on different patients.
- Some patients with very sensitive skin may experience redness under the electrodes after a session. Generally, this redness is completely harmless and usually disappears after 10 to 20 minutes. Never start another stimulation session in the same area, however, if the redness is still visible.
- When using the device to customize programs, take special care that the parameters customized and applied by you to the patient are as you wanted them to be.

**Caution!**

- Do not expose the Cefar URO to direct sunlight, because some of the components may reach unacceptably high temperatures.
- Make sure that pets or vermin do not contaminate the physiotherapy unit and keep them away from it. Also, keep the unit clean and protect it from dust and lint. The safety rules and regulations set forth apply in any case.

**Caution!**

- It is recommended to use the transport bag that comes with the unit, for transport of the device, and to use a proper transport box to ship it.
- Always use the AC adaptor (power supply) provided by the manufacturer to recharge the unit. Never use any other AC Adapter. This could damage the device and generate electrical hazards.
- Only use electrodes supplied by the manufacturer. Other electrodes may have electrical properties that are unsuitable for or may damage the Cefar URO.
- Do not place the electrodes in water.
- Do not apply solvents of any kind to the electrodes.
- **Instructions for electrodes.** See the usage and storage instructions displayed on the bag of electrodes. When the electrodes no longer stick well to your skin, dispose of them in a receptacle out of reach of children and pets
- **Hot casing or battery.** Under extreme use conditions, some parts of the casing might reach up to 118 °F (48 °C). Use caution when manipulating the rechargeable battery right after device use or when holding the device. There is no particular health risk associated with this temperature besides your comfort
- **Size of electrodes.** Do not use electrodes with an active area of less than 10 cm² due to the risk of associated burning. Proceed systematically with caution when the density of the current is over 2 mA/cm².
- If the device is not used for an extended period, we recommend that you charge your device every 5 months. You should store the device in a dry and cool environment.

Note

- For best results, wash and clean the skin of any oil and dry it before attaching the electrodes.
- Never use a set of adhesive electrodes for more than 40 sessions as the quality of the contact between the electrode and the skin, which is essential for the patient's comfort and the effectiveness of the stimulation, gradually deteriorates.
- For information on use and storage please consult the instructions found on the electrodes packaging.

Note**Biocompatibility**

Those parts of the Cefar URO unit that come into contact with the patient when used as intended, are designed to fulfil the biocompatibility requirements of the applicable standards.

4. OVERVIEW

4.1 Device components and accessories

Model: Cefar URO

Part number: 2054

Your kit contains (included in delivery):

Quantity	Description	Part number
1	Cefar URO device	1982
1	Fast charger	2001 (EU), 2133 (UK), 2145 (US)
1	Set of 2 pin cables 2 mm	1990
1	User manual and practical guide	13-01187
1	Transportation pouch	2049
2	Necklace and clip	2023
2	Beltclip	2037

Electrodes are considered as option depending on markets.

4.2 Explanation of symbols (connections and nameplates)



	Follow the instructions of the user manual
	Caution! Observe warnings set forth in operation manual!
	The Cefar URO is a class IIa device with internal electric power and type BF applied parts.
	The power switch On/OFF button is a multi-function button.
	The name and address next to this factory symbol is the manufacturer.

	The number next to this symbol is the article reference number
CE-2797	Device complies with Council Directive 93/42/EEC as amended, concerning medical devices.
	The number next to this symbol is the serial number
	WEEE Mark (European Directive 2002/96/EC). Indicates separate treatment from general waste at end of life.
	Keep dry
IP22	IP classification indicates the degree of protection and thus defines its suitability for use under various ambient conditions.
	Non-ionising radiation
	Not made with natural latex rubber
	Keep away from direct sunlight
	Alternating current input on AC power supply
	Direct current output from power supply
	Protection class II equipment. The AC Power Supply device has double insulation.
	Indoor use only

5. DESCRIPTION OF THE CEFAR URO






5.1 Description of the device

- A ON/OFF
- B Program selection - Pause during stimulation
- C Increase/decrease amplitude increase/decrease timer menu navigation
- D Timer
- E Confirm - Resume stimulation when paused
- F Cables sockets
- G Battery charger socket
- H Belt clip holders
- I Belt clip (can be reversed 180°)
- J Key unlock

5.2 Explanation of symbols

- 1 Battery Status
- 2 Keyboard lock symbol
- 3 Program number (left and right channel)
- 4 Amplitude bargraph (left and right channels)
- 5 Amplitude level in milli Ampere (mA) or Voltage (V)
- 6 Remaining time
- 7 Mixed freq - Muscle stimulation icon
- 8 Pause
- 9 Program end
- 10 Open circuit (disconnect electrode or cable)
- 11 Program Lock

5.3 Battery charge

The Cefar URO is powered by a rechargeable battery.

- A** - To recharge the Cefar URO first disconnect the electrode cables from the device, then plug the charger to a wall socket then connect the stimulator to the charger.
- B** - The battery level of charge is indicated on the treatment screen
- C** - When charging, the battery charge symbol is animated with a battery filling itself from empty to full.
- D** - When battery is fully charge at the end of charging process device will show 100%.

Fig A


Fig B


Fig C - Low charge


Fig D - Half charged


Fig E - Full charge

5.4 Accessories

Electrodes: The Cefar URO should be used with recommended Enovis electrodes

The electrodes will eventually wear out and need to be replaced. It is recommended to replace after approximately 20-40 times of usage.

Note: Stimulation comfort and treatment effectiveness are highly influenced by the quality of the electrodes.

Neck Strap and Clip: The Cefar URO comes with a belt clip and a neck strap/clip allowing you to wear the stimulator around your neck or on your belt and have your hands free during treatment.

To use the neck strap, fix the clip associated to the hole at the top of the unit. (Fig A & Fig B)

The Belt clip can be rotated in order to have the cable outputs on the bottom, or on the top depending on the body area being stimulated. To release the belt clip, push it downwards (Fig C & Fig D).

Fig A

Fig B

Fig C

Fig D

6. STEP BY STEP USE - PATIENT MODE

You must read the contra-indications and safety measures described at the beginning of this manual (Chapter 3 "Safety Information") prior using your stimulator.

6.1 Attach the electrodes


- Connect the electrodes to the cables.
- Attach the electrodes to your body.
- Plug the stimulator cables into the two sockets on the base of the device. The two cables can be connected simultaneously. The cables colours are different to simplify use and facilitate identification of the two channels. Make sure to have the line of the connector upwards to insert the connector into the right direction.
- The kit includes two 2mm pin-connector cables.

6.2 Switch on the stimulator

Press the "On/Off" button.

Note: last program configuration will always be displayed by default.

6.3 Select a preset Program

Press the Program **"P"** button for the left or right channel you want to adjust. The program number is blinking on the screen. Push on up/down arrows of the same channel to increase or decrease the program number and confirm with **"OK"**. If you want to use both channels, then select a program for the other channel. Any combination of programs is possible

Note: When selecting a program the amplitude must be set to 00.0 mA or 0.0 V for the channel concerned. For further information on the programs, see section PRESET PROGRAMS.

6.4 Adjust timer settings if needed

Timer preset value can be either unlimited (---) or with a limited time . See program table summary at the end of this manual. Press once the **"timer"** button, Channel 1 time value is blinking. To adjust time press on adjust value buttons **"up or down"** and press **"OK"** to confirm

Press twice the timer button to select Channel 2 timer adjustment.

6.5 Start the stimulation

Press the **"INCREASE"** button for each channel until you reach a comfortable level of stimulation. Press and hold the button to increase the amplitude gradually.

Note: Always increase the amplitude cautiously!

6.6 Key Lock

An automatically activated key lock prevents unintentional changes during treatment. The key lock is activated if no button is pressed for 10 seconds while stimulation is ongoing. Press any of the **"DECREASE"** buttons to deactivate the key lock.

6.7 Program pause

To pause a program press **"P"**, both Channels 1 and 2 are paused, Amplitude goes down automatically to 00,0. To resume program press **"OK"**. Amplitude increases automatically and progressively to set value.

6.8 Open lead

The broken circuit display indicates on which channel that resistance is too high, or that a cable is broken or detached. A high resistance can be caused by a bad connection between the electrodes and your skin, or that the electrodes need to be replaced. The stimulation is stopped, while trouble shooting is fixed.

Note: Broken circuit is only indicated when running over 9.5 mA. Arrow is pointing to the channel with the open lead.

6.9 Resume treatment after open lead

Press "OK" after checking electrodes and cables connections. Treatment screen is displayed. Amplitudes have to be readjusted manually. A cable breakage can be checked by pressing the cable's pins against one another while increasing the amplitude for the corresponding channel to 11 mA. If the amplitude now drops to 0.0 mA and starts flashing, the cable needs to be replaced.


Note: Never increase the amplitude above 20 mA when you check for cable breaks, since this can damage the stimulator.

6.10 Program end

When the program is finished, the last used programs are selected for the next time the stimulator is turned on.

The stimulator turns off after 5 minutes of inactivity to spare battery.

7. PHYSICIAN MODE

You must read the contra-indications and safety measures described at the beginning of this manual (Chapter 3 "Safety Information") prior using your stimulator.

Physician mode can only be used by trained medical professionals. Other people should not use this mode.

7.1 To access physician mode

When device is switched OFF, press simultaneously "Ch2+" and "P2" buttons and maintain them pressed while pushing on the "ON/OFF" button.

7.2 Physician mode home screen

Navigation using "up and down", press "OK" to confirm

Lock Program: The Cefar URO can be locked to prevent accidental changing of programs

Custom Program: With the Cefar URO it is possible to create and store three custom programs for individual needs

Access Compliance: Compliance gives you the possibility to monitor the use of the stimulator

Reset Compliance: Reset compliance data

Set CC-CV: Cefar URO can operate on Constant Voltage or Constant current mode

Enter Pain score: Patient pain can be assessed with numerical scale

Reset pain score: Patient pain score can be reset

7.3 Lock Program

It is possible to lock a program to avoid using an incorrect program. After selecting a program in patient mode, use up and down to select **"Lock/Unlock Program"**, **"OK"** to confirm and go back to Physician home screen.

7.4 Custom Programs


7.4.1 Select Custom Programs

CP1, CP2, CP3 can be customized for individual needs.

Use **"up and down"** to select program, **"OK"** to confirm.

7.4.2 Custom Program first screen

Stimulation mode can be set as Conventional, Burst, Pulse Width Modulation (PWM), Frequency Modulation, Mixed Frequency. Press **"up and down"** to select appropriate mode in the parameters range. Press **"OK"** to confirm.


7.4.3 Summary screen

Summary screen is displayed after last parameter is adjusted and "OK" button is pressed once. Use "up and down" to scroll through the different parameters for visualization only.

PC2	
Conventional	
Frequency	80
Pulse duration	180

7.4.4 Parameters range

7.5 Compliance

Device displays session time as time of use since last reset. Total time of device use can't be reset. Session time can be reset.

Session time **2h35m**
Total time **142h 19m**

Compliance

Reset

Cancel

7.6 CC-CV

Cefar URO can be set on Constant Current or Constant Voltage mode Press "up and down" to select appropriate mode. "OK" to confirm and back in Physician home screen


CC - Constant current: This mode guarantees optimal stimulation without variation of amplitude by automatically adapting the voltage when resistance varies during the stimulation. The resistance may vary depending on tissue impedance, or electrodes contact quality. The CC mode is recommended by default.

CV - Constant Voltage available only on 1 channel: This mode adapts automatically the amplitude when the resistance varies, allowing more comfortable stimulation if the patient loses electrode surface contact or takes electrodes the electrodes in their fingers while stimulation is on.

CV mode is recommended to be used when moving electrodes from place to place or when the patient does not feel comfortable managing the device.

Constant Current

Constant Voltage

7.7 Pain score

- Pain score can be assessed and stored in a numerical scale from "0" absence of pain to "10" maximum pain.
- Use "Up and Down" to set pain score
- Press "OK" to confirm
- Pain score can be reset

Pain score (0-10)

5

Pain score **5**

Reset

Cancel

8. PROGRAMS

8.1 Program parameters

Program Number	Name	Pulse Width 1	Frequency 1
P1	NTP 15Hz	150 µs	15 Hz
P2	NTP 14Hz	210 µs	14 Hz
P3	NTP 10Hz	200 µs	10 Hz
P4	NTP 50Hz	200 µs	50 Hz
CP1	Custom program 1	-	-
CP2	Custom program 2	-	-
CP3	Custom program 3	-	-

8.2 Program description and recommendations

The choice of program is made according to the prescriber's recommendation.

Depending on the chronaxis point and the sensitivity of the patients, the P1 to P4 programs offer configurations that provide satisfaction.


The electrodes are placed on the path of the posterior Tibial Nerve. On the medial side (inside) of the ankle, we recommend placing one electrode in the submalleolar area and another in the retro-malleolar area. Be careful, the electrodes must never touch each other. Set the amplitude so that the stimulation produces intense, painless tingling (paresthesia) in the arch of the foot.

The CP1, CP2, and CP3 programs allow prescribers to customize the stimulation parameters to best suit the patient's needs.

9. TROUBLESHOOTING

9.1 Electrode or cable failure

The Cefar URO will display a pair of electrodes and an arrow pointing to the channel where the problem was detected.

The electrode error message can mean that:

- There are no electrodes connected to this channel.
- The electrodes are old, worn, and/or the contact is poor: try using new electrodes.
- The electrode cable is faulty: try connecting it to another channel. If the problem persists, replace the cable.

If for any reason your device seems to present a malfunction other than those presented above, please contact Enovis customer service.

Problem	Possible Cause	Solution
Electrode failure	Poor electrode connection to the cable	Check that the electrodes are properly connected to the cable.
	Poor electrode connection to the skin	Check if the electrodes are outdated, worn and/or the contact is poor: try with new electrodes

9.2 Battery Level

Problem	Possible Cause	Solution
Device battery low	During the stimulation the device may turn out due to discharged battery.	Stop the stimulation and recharge the device.

9.3 Others

Problem	Possible Cause	Solution
Display does not come on	Low battery	Charge battery
Weak stimulation with loaded batteries	Electrodes dried out, lost their adhesive power and have no adequate connection to the skin.	Replace electrode
	Electrode placement	Make the electrodes at least 5 cm (2 inches) apart
Stimulation stops with loaded batteries	Poor electrode contact	Reapply electrodes, secure firmly. Electrodes must be a minimum of 5 cm apart.
	Damaged or worn electrodes	Replace
	Cabel or contacts damaged	Replace
Stimulation weakens within minutes of starting treatment with loaded batteries	This is a normal body adaptive process	Increase the amplitude (intensity) if required.
Stimulation is uncomfortable	Amplitude (intensity) is too high	Decrease amplitude (intensity).
	Electrodes are too close together	Reposition the electrodes. Electrodes must be a minimum of 5 cm (2 inches) apart.
	Damaged or worn electrodes	Replace
	Ensure proper program is being used	A. Refer to the programs description B. Contact clinician if discomfort persists
Stimulation is ineffective	Improper electrode placement	Reposition electrodes. Electrodes must be a minimum of 5 cm (2 inches) apart.
	Unknown	Contact clinician.

Stimulation only felt on one electrode	Improper electrode placement	A. Reposition electrode. Electrodes must be a minimum of 5 cm (2 inches) apart. B. Replace electrodes.
Stimulation on one channel (side) only	Electrodes A. Worn or damaged B. Improper placement	A. Replace B. Reposition electrode. Electrodes must be a minimum of 5 cm (2 inches) apart.
Intermittent Output	Intermittent Program in use	Some programs will seem intermittent. This is expected. Refer to section 8.1 for a description of the programs.
Stimulation is not producing the usual sensation	Settings and Electrodes positioning	A. Check that all the settings are correct and ensure the electrodes are positioned properly B. Change the positioning of the electrodes slightly.

When an error is detected by the device (from anywhere in the device), the stimulation shall stop on all channels and a specific screen will be displayed. If this screen is displayed, record the Error number and call the Enovis Customer Service.

10. CARE, MAINTENANCE, TRANSPORT, ENVIRONMENTAL STATEMENT

Warning!

Shock hazard - Remove the power cord of the device from the wall outlet before cleaning.
Shock hazard, equipment damage -

- Liquids must not enter the device and its components. If liquids have entered into the components, the device must be immediately checked by a service technician, before it can be reused.
- Never dismantle the device or the AC adapter as they contain high-voltage parts with a risk of electric shock.

- All parts of the device can be disinfected weekly by wiping down with a disinfectant. Thus, it complies with the special hygiene standards for medical technical equipment.
- All components can be cleaned with common disinfectants and mild household detergents.
- Only use a soft cloth and an alcohol-based, solvent-free cleaning product, to wipe the therapy unit down.
- Allow the device to completely dry before use.
- Please refer to environmental part (11.1.4) for storage and transport conditions.

Warning!

Patient hazard - patient contamination

- Before using the unit on another patient, clean and disinfect it according to the instructions in this section.

Caution!

Equipment damage -

- The plastic material used is not resistant to mineral acids, formic acid, phenols, cresols, oxidants and strong organic or inorganic acids with a pH value below 4.
- Use only clear disinfectants to prevent discolouration of the device.
- Do not expose the therapy unit to strong ultraviolet radiation (sunlight) and fire.
- Do not sterilize the stimulator.
- Do not immerse in liquids.

10.1 Maintenance

Your device does not require calibration or frequently safety testings. Each stimulator is tested prior to distribution. Its characteristics do not vary under normal conditions. If your stimulator contains parts that seem worn or defective, please stop using it and contact the customer service centre that has been stipulated and authorised by the manufacturer regarding an upgrade.

There are no user serviceable parts inside the device. Please contact Enovis if you need assistance setting up, using or maintaining the Cefar URO or to report any unexpected operations or events. Please also contact Enovis for clinical questions or if any component of the Cefar URO is not functioning properly or requires servicing. You can find the different contact possibilities under enovis-medtech.eu

Warning!

Shock hazard, Equipment damage -

Do not attempt to repair the stimulator or any of its accessories. Never dismantle the device because of risk of electric shock. Enovis declines all responsibilities for any damages or consequences resulting from unauthorised attempts to open, modify, or repair the stimulator. This may only be done by persons or repair services authorised by the manufacturer.

10.2 Transport

10.2.1 Transport of the Cefar URO

1. Prepare the device and its accessories for shipping within the original Cefar URO pouch
2. Turn the device and it's accessories off.
3. Disconnect and dismount the device and it's accessories by following the guidelines
4. Place the accessories within the box
5. Store the user manual in the pocket of the transport bag.

Caution!

Equipment damage -

Only use the original transport bag for carrying the device around.

**Caution!**

Equipment damage -

Only use the original shipping box for shipping the device. Enovis cannot be held liable for transport damage if the device is not packed in its original shipping box.

10.3 Environmental statement, expected life

The device is electronic equipment and may include substances that can damage the environment. It must not be disposed of with unsorted household or municipal waste. It requires separate disposal at a suitable collection point for recycling of electronic equipment. By doing so, you will be contributing to the safeguarding of natural resources and health. Please contact Enovis for information about the possible recycling of the product

Refer to the TRIMAN logo applicable in France for sorting appliances, accessories and packaging.

The product as well as the parts and accessories supplied with it are designed for a minimum service life of 5 years of normal usage.

11. TECHNICAL DATAS, STANDARDS, GUARANTEE, PATIENTS

11.1. Technical data

11.1.1 General information

Model: Cefar URO

Part number: 1982

EU MDR: Class IIa

IP class: IP22

Applied part: Type BF

Battery: Lithium Polymer (Li-Po) rechargeable battery (3.7V 1500 mAh).

Battery chargers: Only battery chargers provided by Enovis™ and with P/N 2001 (EU), 2133 (UK), 2145 (US) can be used to recharge the batteries supplied with Cefar URO stimulators.

Weight: 156g

Dimensions: 132 x 60 x 24 mm

11.1.2 Stimulation parameters

All electrical specifications are given for an impedance of 1500 ohms per channel (60mA) and for 1000 Ohms per channel for 99.5 mA.

Outputs: Two independent and individually adjustable channels that are electrically isolated from each other.

Pulse shape: Constant asymmetric rectangular current with pulse compensation for TENS programs to eliminate any direct current component to prevent residual polarisation at skin level.

Maximum pulse intensity: 99.5 mA in Constant Current mode or 99 V in Constant Voltage mode for neurostimulation.

Pulse intensity increments: Manual adjustment of stimulation intensity from 0 to 99 mA in minimum increments of 0.5 mA for Constant Current Mode, 0 to 99 V in 1 V increments for Constant Voltage mode for electrostimulation.

Pulse width: 60 to 400 µs.

Frequency: 1-120 Hz

Maximum electrical charge per pulse: 80 micro coulombs (2 × 40 µC, compensated)

11.1.3 Information on electromagnetic compatibility (EMC)

The Cefar URO is designed to be used in typical environments that have been approved in accordance with the EMC safety standard EN 60601-1-2.

This device complies with the CISPR standard, indicating that radio frequency (RF) emissions are not likely to cause interference with electronic equipment installed nearby (radios, computers, telephones, etc.).

The Cefar URO is designed to withstand foreseeable disturbances from electrostatic discharge, magnetic fields from the mains power supply or RF transmitters.

Nevertheless, it is not possible to ensure that the stimulator will not be affected by mobile powerful RF (radio frequency) fields from other sources.

For more detailed information concerning electromagnetic emissions and immunity, refer to the EMC tables.

11.1.4 Environmental conditions

Storage and Transport Conditions

The device must be stored and transported in accordance with the following conditions:

Temperature: -20°C to 70°C

Maximum relative humidity: 75%

Atmospheric pressure: 700 hPa to 1,060 hPa

Conditions of use

Temperature: 10°C to 40°C

Maximum relative humidity: 30% to 75%

Atmospheric pressure: 700 hPa to 1,060 hPa

11.2 Standards

The Cefar URO has been designed, manufactured, and distributed in compliance with the requirements of EU MDR, as amended, on medical devices.

Current international standards require that a warning be given concerning the application of \ electrodes to the thorax (increased risk of cardiac fibrillation).

The Cefar URO also complies with Directive 2012/19/EU on waste electrical and electronic equipment (WEEE).

11.3 Guarantee

This guarantee is valid only if it is accompanied by proof of purchase.

Your statutory rights are not affected by this guarantee.

Your Cefar URO stimulator is guaranteed for a period of 3 years from the date of purchase. The guarantee covers the device and the AC adaptor (hardware and labour), but not the cables, battery and electrodes.

All defects resulting from poor quality material or workmanship are covered. This guarantee does not cover damage resulting from impact, accidents, misuse, inadequate protection against moisture, immersion in water or repairs made by unauthorised personnel.

12. FREQUENTLY ASKED QUESTIONS

FOR HOW LONG TIME CAN I STIMULATE?

Follow the recommendations given by your health care provider. In general, the posterior tibial nerve TENS method can be performed without time limit, but a stimulation session should last at least 20 minutes.

FOR HOW LONG LAST THE ELECTRODES?

The self-adhesive electrodes last for approximately 20 to 40 occasions. The durability depends on how good the care and maintenance instructions are followed. Read the care and maintenance instructions delivered with the electrodes.


13. CONTACTS

We would be happy to answer any questions you may have about our products and services. Please contact your local dealer, or your corresponding Enovis site. Enovis sites are listed on the backside of the cover.

For technical service from Enovis, please contact:

internationalproductsupport@enovis.com

14. ELECTRODE PLACEMENT GUIDE

Use two round electrodes, 32 mm diameter, on a single internal ankle and place them as :

- the first retro-malleolar
- the second in the submalleolar

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5191 Sea Otter Place, Suite 200, Carlsbad, CA 92010, U.S.A

Chattanooga Rehab.fr

Individual results may vary. Neither Enovis, DJO, LLC or any of their subsidiaries dispense medical advice. The contents of this document do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.

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DJO France SAS