

# MEDICAL OXYGEN CONCENTRATOR

MODEL: JAY-10



**DUAL FLOW**  
USA imported humidifier  
6PSI pressure release



Innovative “cooling”  
technology to improve  
reliability and life of the unit



**Big LCD display**  
24H\*365D work continuously

## Optional

1. Nebulizer(Atomization)
2. Low purity alarm with low flow alarm
3. High temperature alarm
4. Pulse oximeter
5. Dual flow

## Specification

Flow Rate	0-10LPM
Purity	93% ( $\pm 3\%$ )
Outlet Pressure	0.04-0.07MPA (6-10PSI)
Sound Level	$\leq 50$ db
Power Consumption	$\leq 880$ w
LCD Display	Switch times, operating pressure, present working time, accumulated working time, presetting time from 10mins to 40hours
Alarm	Power failure alarm, High&low pressure alarm, Low purity alarm
Net Weight	27 Kgs
Size	365×375×600mm